



## STEM CELL LABORATORY (STCL)



**DOCUMENT NUMBER:** FLOW-GEN-049 JA4

**DOCUMENT TITLE:**

Color Compensation Settings for BD FACSCanto II

**DOCUMENT NOTES:**

### Document Information

**Revision:** 01

**Vault:** FLOW-General-rel

**Status:** Release

**Document Type:** FLOW

### Date Information

**Creation Date:** 06 Jun 2025

**Release Date:** 13 Jun 2025

**Effective Date:** 13 Jun 2025

**Expiration Date:**

### Control Information

**Author:** WATER002

**Owner:** WATER002

**Previous Number:** None

**Change Number:** STCL-CCR-565

## FLOW-GEN-049 JA4

### Color Compensation Settings for BD FACSCanto II

Follow this procedure to reset compensation settings after service/maintenance, if compensation values are out of range, or when checking daily QC if populations deviate from original gates.

**NOTE:** Will need 12 x 75 mm plain tubes (***NOT*** *TruCount tubes*)

- **Tube C-1:** Unstained
- **Tube C-2:** Add 10uL CD3 FITC (BD Cat. No. 349201)
- **Tube C-3:** Add 2.5uL CD3 PE (BD Cat. No. 340662)
- **Tube C-4:** Add 7uL CD3 PerCP (BD Cat. No. 340663)
- **Tube C-5:** Add 2.5uL CD8 APC (BD Cat. No. 340659)
- **Tube C-6:** Add 15uL 7AAD only (BD Cat. No. 559925 or 51-68981E)

#### **Staining Protocol:**

- Add appropriate amounts of antibody to each respective tube (*amounts included above*)
- From Status Flow Cells, Lymphs
- Add 50uL to each tube
- Vortex samples
- Incubate at room temperature (*in the dark*) for 15 minutes
- Lyse, using 450uL/tube, with PharmLyse (BD Cat. No. 555899; using 10% stock solution
  - To prepare the 10% stock solution, add 1-part PharmLyse to 9-parts of sterile, distilled, filtered water from the deionized water system labeled “**spigot # 1**”.
- Incubate at room temperature (*in the dark*) for 15 minutes

#### **Acquisition Protocol:**

##### **Diva software:**

- Acquire in Comp 4-Color Template for 7AAD or PerCP (*these fluorochromes are in same channel*)
- Apply application settings before running comps
- Link to previous comp experiment
- Must acquire using Global Worksheet (WS)
- Normal WS will automatically be created when “create compensation controls”
- When start new CST bead lot, we must update application settings, which could affect comps
- Check the cytometer settings by adding an extra specimen tube (not part of comp controls) and don’t record, just acquire a few events for each parameter to see if population looks good and is on scale
- Check FSC/SSC settings at beginning of run before setting up carousel; otherwise, will run whole carousel at settings for first tube before they’re adjusted
- Save new comp settings and overwrite the old 4-color comp template with new ones (**NOTE:** It is not recommended by BD to try and save the old settings as this could cause confusion in choosing the incorrect settings)
- Check settings daily when acquiring Status Flow Controls.

**Signature Manifest****Document Number:** FLOW-GEN-049 JA4**Revision:** 01**Title:** Color Compensation Settings for BD FACSCanto II**Effective Date:** 13 Jun 2025

All dates and times are in Eastern Time.

**FLOW-GEN-049 JA4 Color Compensation Settings for BD FACSCanto II****Author**

Name/Signature	Title	Date	Meaning/Reason
Barbara Waters-Pick (WATER002)		11 Jun 2025, 06:25:20 PM	Approved

**Management**

Name/Signature	Title	Date	Meaning/Reason
Barbara Waters-Pick (WATER002)		11 Jun 2025, 06:25:34 PM	Approved

**Medical Director**

Name/Signature	Title	Date	Meaning/Reason
Joanne Kurtzberg (KURTZ001)		11 Jun 2025, 06:32:37 PM	Approved

**Quality**

Name/Signature	Title	Date	Meaning/Reason
Melissa Ritt (MSR68)	GMP, Quality Assurance Associate I	12 Jun 2025, 09:16:52 AM	Approved

**Document Release**

Name/Signature	Title	Date	Meaning/Reason
Amy McKoy (ACM93)	Document Control Specialist	12 Jun 2025, 09:29:00 AM	Approved