

STEM CELL LABORATORY (STCL)



DOCUMENT NUMBER: STCL-PROC-022 FRM1

DOCUMENT TITLE:

Stem Cell Laboratory, STCL, Clinical Hematopoietic Progenitor Cell Assay (HPCA) Worksheet

DOCUMENT NOTES:

Was STCL-SOP-029

Document Information

Revision: 03 Vault: STCL-Processing-rel

Status: Release Document Type: Processing

Date Information

Creation Date: 27 Jan 2014 Release Date: 10 Mar 2014

Effective Date: 10 Mar 2014 Expiration Date:

Control Information

Author: WATE02 Owner: WATE02

Previous Number: STCL-PROC-022 FRM1 Rev | Change Number: STCL-CCR-168

STCL-PROC-022 (FRM 1) Stem Cell Laboratory (STCL) Clinical Hematopoietic Progenitor Cell Assay (HPCA) Worksheet

Recipient Name:					Recipient History #:					
Donor name:					Donor History #:					
Date/Time received:					Date/Time Collected:					
(Check	one): ☐Fr	esh F	IPCA = 2 x10 ⁵ ce	ells/ml (Pla	ting 1 x 10 ⁴ co	ells/well)	☐ Frozen HI	$PCA = 4 \times 10^5$	cells (Plating 2	2 x 10 ⁴ cells/well)
Date	Source (CBU,PB,P BM, Othe	SC,	Test Requested (HPB, HPA, etc		cession ‡	#	Cell Count x 10 ⁵ /ml	Volume needed for test	Tech Initials	Verified by (Tech Initials)
CALCUI	ATIONS						1978			15.9
Sou	Source Raw Cell Count (CBU,PB,PSC,BM,Other) x 10 ⁶ /ml		Via			tal Specimen ume Counted			Tech Initials	
						-				

			etal. 7.0				-			
RESULT	rs		1010							
Source		M	GEMM	BFUE		d by: nitials)	Date			
						15				
				8. ***						
				120			= = = = =			
Comments	:									·
					2					
	1,000									

INSTRUCTIONS

Field	Requirements		
Recipient Name	Enter recipient's name		
Recipient History #	Enter recipient's medical history number of patient/donor		
Donor Name	Enter donor's name		
Donor History #	Enter donor history #		
Date/Time received	Enter date and time unit received into laboratory		
Date/Time collected	Enter date and time unit was collected		
Check (ONE) Fresh HPCA or Frozen HPCA	Fresh products require 2 x 10e5 cells/ml (in order to plate 1x10e4 cells/well; Frozen product require 4 x 10e5 cells/ml (in order to plate 2x10e4 cells/well).		
Date	Enter the date the specimen was plated.		
Source Enter cell source being tested (i.e. cord blood (CBM), peripheral blood (PB). etc.			
Test	Enter test being performed on unit (ie. HPB, HPA)		
Accession #	Enter accession number generated by Pathnet (when ordered)		
Cell count/ml	Cell count of product (x 10e5 cells per ml)		
Volume needed for test	Volume of sample needed to plate 2x10 ⁵ cells or 4 x 10e ⁵ depending on whether the sample is fresh or thawed.		
Tech Initials	Initials of technologist entering information on worksheet		
Verified by (Tech Initial)s	Initials of lab tech verifying the calculations and the initials of the		
~	tech entering the results in Pathnet		
Source	Enter type of specimen collected (i.e. cord blood unit (CBU),		
Daniel and	bone marrow (BM), peripheral blood (PB)		
Raw cell count	Cell count of product (x 10e6 cells / ml)		
Viability	Enter percent viability by trypan blue (# of live cells)		
Total Specimen Volume counted	Total volume of sample counted		
Volume of sample used QS to 1.0 ml	Volume used to set up HPC assay		
Tech Initials	Initials of tech performing calculations		
Results	Record the sample source and the number of GM, GEMM and BFUE colonies there were counted on the plate. Results for each colony type will be added together and an average obtained. Enter the tech initials and date the plate was read / counted.		
	RESULTS: (Fresh) – Multiply the average colony count by 10 and then enter the result for each colony type in Pathnet. (Frozen / Thawed) – Multiply the average colony count by 10 and then divide by 2; Enter the result for each colony type in Pathnet.		
Comments	Enter any comments associated with this clinical specimen.		

Signature Manifest

Document Number: STCL-PROC-022 FRM1

Revision: 03

Title: Stem Cell Laboratory, STCL, Clinical Hematopoietic Progenitor Cell Assay (HPCA)

Worksheet

All dates and times are in Eastern Time.

STCL-PROC-022 FRM1 Stem Cell Laboratory, STCL, Clinical Hematopoietic Progenitor Cell Assay (HPCA)

Author				
Name/Signature	Title	Date	Meaning/Reason	
Barbara Waters-Pick (WATE02)		31 Jan 2014, 07:14:12 PM	Approved	
Manager				
Name/Signature	Title	Date	Meaning/Reason	
Barbara Waters-Pick (WATE02)		31 Jan 2014, 07:14:25 PM	Approved	

	D : .
MACHICAL	Director

Name/Signature	Title	Date	Meaning/Reason
Joanne Kurtzberg (KURTZ001)		03 Feb 2014, 09:06:45 PM	Approved

Quality

Name/Signature	Title	Date	Meaning/Reason
John Carpenter (JPC27)		04 Feb 2014, 11:15:39 AM	Approved

Document Release

Name/Signature	Title	Date	Meaning/Reason	
Sandy Mulligan (MULLI0	26)	28 Feb 2014, 05:59:50 PM	Approved	