



STEM CELL LABORATORY (STCL)



DOCUMENT NUMBER: STCL-EQUIP-005 JA1

DOCUMENT TITLE:

Review of Control Rate Freezing Graphs JA1

DOCUMENT NOTES:

Document Information

Revision: 01

Vault: STCL-Equipment-rel

Status: Release

Document Type: Equipment

Date Information

Creation Date: 21 May 2014

Release Date: 23 Jun 2014

Effective Date: 23 Jun 2014

Expiration Date:


Control Information

Author: WATE02

Owner: WATE02

Previous Number: None

Change Number: STCL-CCR-133




Step 1 – When selecting the RUN button the FIRST time, the chamber of the CRF will maintain a temperature of ~ 4° until RUN has been pressed a SECOND time. This pre-cools the chamber to be cooled before bags are placed in the CRF to start the freezing run.


Step 2 - Once RUN has been pressed the SECOND time, the temperature will drop 25°C/min until the **chamber** reaches -40°C.

Step 3 – Temperature will drop 50°C/min until the **chamber** reaches 0°C.

Step 4 – **HOLD Chamber** at 0°C for 15 minutes.




Step 5 – Temperature will drop 1°C/min until **sample** reaches -10°C.




Step 6 – Temperature will drop 30°C/min until **chamber** reaches -65°C

Step 7 – Temperature will drop 15°C/min until **chamber** reaches -20°C.

Step 8 – Temperature will drop 1°C/min until **chamber** reaches -70°C (this step takes approximately 50 minutes given the temperature drop of 1°C/min.




Step 8 – Continuation of step 8 (each horizontal line represents 1 minute (this screen shot spans a total of approximately 26 minutes)).



Step 9 – Temperature will drop 10°C/min until **chamber** reaches -90°C.

Step 10 – **HOLD Chamber** for 6 minutes at -90°C



END - Holds **chamber** at **-90°C** until the CRF is stopped. Chamber stays at **-90°C** until **BACK** is pressed. Do **NOT** press **BACK** until you are ready to remove cryopreserved products from the Control Rate Freezer for storage in designated LN2 freezer.

TROUBLESHOOTING:

- If no Heat of Fusion (HOF) is reflected on the freezing graph, a segment from one bag from each recipient whose product was included in that freezing run must be tested for viability. If thawed viability is <85%, consult with the attending physician in case there they want to collect additional cells from the donor/recipient. (**NOTE:** Sometimes when there is a low TNC or the temperature probe is not placed properly between two bags, this can result in HOF not being displayed properly on the graph).
- File the reviewed freezing graphs indefinitely in the designated file cabinet in the laboratory and/or in Record Retention department once on-site storage space has been exhausted.

Signature Manifest**Document Number:** STCL-EQUIP-005 JA1**Revision:** 01**Title:** Review of Control Rate Freezing Graphs JA1

All dates and times are in Eastern Time.

STCL-EQUIP-005 JA1 Review of Control Rate Freezing Graphs**Author**

Name/Signature	Title	Date	Meaning/Reason
Barbara Waters-Pick (WATE02)		22 May 2014, 07:53:17 AM	Approved

Manager

Name/Signature	Title	Date	Meaning/Reason
Barbara Waters-Pick (WATE02)		22 May 2014, 07:53:46 AM	Approved

Medical Director

Name/Signature	Title	Date	Meaning/Reason
Joanne Kurtzberg (KURTZ001)		24 May 2014, 06:33:24 AM	Approved

Quality

Name/Signature	Title	Date	Meaning/Reason
John Carpenter (JPC27)		27 May 2014, 12:31:47 PM	Approved

Document Release

Name/Signature	Title	Date	Meaning/Reason
Sandy Mulligan (MULLI026)		11 Jun 2014, 03:40:53 PM	Approved