



## STEM CELL LABORATORY (STCL)



<b>DOCUMENT NUMBER:</b> STCL-PROC-035 FRM1
<b>DOCUMENT TITLE:</b> Manual Differential, Cord Blood Units
<b>DOCUMENT NOTES:</b>

### Document Information

<b>Revision:</b> 05	<b>Vault:</b> STCL-Processing-rel
<b>Status:</b> Release	<b>Document Type:</b> Processing

### Date Information

<b>Creation Date:</b> 22 May 2012	<b>Release Date:</b> 04 Jun 2012
<b>Effective Date:</b> 04 Jun 2012	<b>Expiration Date:</b>

### Control Information

<b>Author:</b> WATE02	<b>Owner:</b> WATE02
<b>Previous Number:</b> 9D.211.03 (FRM1) STCL-PRC <b>Change Number:</b> STCL-CCR-034	

**STCL-PROC-035 (FRM1)**  
**Manual Differential (Slide Method) - CBUs for CCBB Program**

Manual 100-Cell Differential (Cord Blood Units for CCBB Program) – (%) Percentage

<b>Bar Code</b>
-----------------

**Blasts****NRBC****Promyelocytes****Myelocytes****Metamyelocytes****Bands****Segmented Neutrophils****Eosinophils****Basophils****Lymphocytes****Monocytes****Plasma cells****Other**Platelet Est / Morph? ☐ Yes ☐ No

<b>PRE-PROCESSING</b> % Collection bag (Pre-Hespan)	<b>POST-PROCESSING</b> % (Volume depleted/RBC reduced)
Plt Estimate _____ K	Plt Estimate _____ K
Morphology OK? ____ Y ____ N	Morphology OK? ____ Y ____ N

% NRBCs \_\_\_\_\_ / \_\_\_\_\_      % Mononuclear cells \* \_\_\_\_\_ / \_\_\_\_\_

\* Mononuclear Cells = Blasts + Lymphs + Monos

\_\_\_\_\_  
**Technologist's signature**

\_\_\_\_\_  
**Date**

\_\_\_\_\_  
**Study ID**

*N/A = Not Applicable*

STCL-PROC-035 (FRM1) Manual Differential (Slide Method) - CBUs for CCBB Program

Example

STCL, DUMC

Durham, NC

INSTRUCTIONS	
Field	Requirement
<b>PRE-PROCESSING</b> % Collection bag (Pre-Hespan)	Enter the raw count of each cell type and their percentage for the collection bag (when counting 100 cells). For CCBB UCB units, the NRBC count is INCLUDE NRBCs in the 100 cell differential.
<b>POST-PROCESSING</b> % (Volume depleted/RBC reduced)	Enter the raw count of each cell type and their percentage for the collection bag (when counting 100 cells). For CCBB UCB units, the NRBC count is INCLUDE NRBCs in the 100 cell differential.
<b>Platelet Est / Morph?</b> <input type="checkbox"/> Yes <input type="checkbox"/> No	<p>If <b><u>YES</u></b> is checked, perform platelet estimate and morphology on the slide.</p> <ul style="list-style-type: none"> <li>• Platelet Estimate should be entered on the line provided (ie. 20,000 = 20K)</li> <li>• Morphology OK? (Check) YES if morphology is OK, (Check) NO, if morphology is NOT normal.</li> </ul> <p>If <b><u>NO</u></b> is checked, leave fields blank as no platelet estimate or morphology check is being requested.</p>
<b>% NRBCs</b>	Enter the raw # of NRBCs from the Pre-Processing fraction and the Post-Processing fraction (Pre # / Post #).
<b>% Mononuclear Cells (MNCs)</b>	Enter the raw # of MNCs from the Pre-Processing fraction and the Post-Processing fraction. The % MNCs is equal to the # of Blasts + Lymphs + Monos.

STCL-PROC-035 (FRM1) Manual Differential (Slide Method) - CBUs for CCBB Program

Example

STCL, DUMC

Durham, NC

## Signature Manifest

Document Number: STCL-PROC-035 FRM1

Revision: 05

Title: Manual Differential, Cord Blood Units

All dates and times are in Eastern Time.

## STCL-PROC-035 FRM1 Man Diff CCBB

## Author Approval

Name/Signature	Title	Date	Meaning/Reason
Barbara Waters-Pick (WATE02)		22 May 2012, 05:46:43 PM	Approved

## Manager Approval

Name/Signature	Title	Date	Meaning/Reason
Barbara Waters-Pick (WATE02)		22 May 2012, 05:47:03 PM	Approved

## Medical Director Approval

Name/Signature	Title	Date	Meaning/Reason
Joanne Kurtzberg (KURTZ001)		22 May 2012, 07:05:16 PM	Approved

## QA Approval

Name/Signature	Title	Date	Meaning/Reason
Linda Sledge (SLEDG006)		23 May 2012, 08:31:38 AM	Approved

## Notification

Name/Signature	Title	Date	Meaning/Reason
Barbara Waters-Pick (WATE02)		23 May 2012, 08:31:38 AM	Email Sent

## Review: STCL-PROC-035 FRM1 05

## Review

Name/Signature	Title	Date	Meaning/Reason
Barbara Waters-Pick (WATE02)		05 May 2014, 07:20:24 AM	Reviewed
Joanne Kurtzberg (KURTZ001)		05 May 2014, 10:12:33 AM	Reviewed
(DOHNA02W) for John Carpenter (JPC27)		05 May 2014, 12:23:51 PM	Reviewed
Betsy Jordan (BJ42)		06 May 2014, 10:35:03 AM	Reviewed
Sharon Hartis (SH259)		06 May 2014, 10:39:22 AM	Reviewed