



STEM CELL LABORATORY (STCL)



DOCUMENT NUMBER: STCL-PROC-037 FRM1

DOCUMENT TITLE:

Miltenyi Large Scale CCB \(\text{IFN-gamma}\) Worksheet FRM1

DOCUMENT NOTES:

Document Information

Revision: 02

Vault: STCL-Processing-rel

Status: Release

Document Type: STCL-Processing

Date Information

Creation Date: 11 Jul 2012

Release Date: 13 Dec 2012

Effective Date: 13 Dec 2012

Expiration Date:

Control Information

Author: WATE02

Owner: WATE02

Previous Number: STCL-PROC-037 FRM1 Rev 1 **Change Number:** STCL-CCR-048

Barcode

**STCL-PROC-037 FRM1
MILTENYI LARGE SCALE CCS (IFN-gamma) WORKSHEET**

Donor Name: _____ History No.: _____
 Recipient Name: _____ History No.: _____
 Collection Date and time: _____ Bag labels checked by: _____
 Selection Date: _____ Performing technologist: _____

Pre Selection

Centrifuge turned off by: _____ Time: _____ RPMI brought to room temperature by: _____
 AB Serum thawed by: _____ AB Serum filtered by: _____
 RPMI with 2% by volume filtered AB serum prepared by: _____
 mls RPMI: _____ mls AB serum: _____ Transferred to 600ml transfer pack placed in
 refrigerator by: _____
 (600mls will be needed)
 RPMI with 10% by volume filtered AB serum prepared by: _____
 mls RPMI: _____ mls AB serum: _____ Transferred to 600ml transfer packs and placed
 in 37° incubator by: _____
 (1200mls will be needed)
 R.T. RPMI transferred to 1000 ml bag: _____ Transferred by: _____
 (Transfer by aseptically pouring RPMI through 60cc syringe attached to luer adapter)
 3 liters of working buffer prepared, labeled and placed in the refrigerator by: _____
 300ml transfer pack prepared and weighed for final cell collection by: _____
 600ml transfer packs labeled appropriately by: _____
 Hood inspected, cleaned and disinfected prior to initiating procedure by: _____

Initial Product (Sample A)

Cell count: _____ Volume: _____ Total Cells: _____ Viability: _____
 Product Cultured: _____ Flow sent _____ HPCA sent _____ ABO: _____ by _____
 Calculations for 1×10^9 cells: _____

Cell Processing - Washing

Inject 1×10^9 cells into 600 ml transfer pack through a luer adapter and fill bag to 500cc with R.T.
 RPMI. Centrifuge at R.T. at 200xg for 10 minutes with no brake. Express supernatant.
 Perform a second wash of the 1×10^9 cells in the same manner.
 First wash completed by: _____ Second wash completed by: _____
 Post 2nd wash, supernatant expressed to minimum volume by _____ Volume = _____
 (target=40mls)
 Volume of 37°C RPMI with 10% AB serum to add to cell prep:
 (100mls – volume remaining in bag) = _____ Calculated and added by _____

Barcode

Cell Processing – Re-stimulation

Transfer to a gas-permeable cell culture bag and add 1 ml of pp65. Added by _____

Incubate for 4 hours at 37°C. Time incubation started: _____ Timer set: _____

Note: Target is 100ml of the product with a cell density of 0.5 to 1.0×10^7 cells/cm² and 1×10^7 cells/mL.**During Incubation**

Fill three 600 ml transfer packs with 125 mls warm RPMI with 10% AB serum, label as #1, #2 and #3 and return to incubator. Filled, labeled and stored by: _____

Verify a chilled (4°C) centrifuge is ready.**Cell Processing – Labeling (IFN-gamma Catchmatrix Reagent)**

Product removed from incubator by: _____ Product transferred to a 600ml transfer pack and filled with cold (2-8°C) RPMI with 2% AB serum added by: _____

Product centrifuged at 4-8°C at 300xg for 10 minutes with no brake by: _____

(During centrifugation, prepare 4 buckets of ice).

Supernatant expressed and cells re-suspended by: _____

Volume of cells remaining: _____ (target is smallest volume feasible)

Catch Matrix injected into cell mixture by: _____ (20g needle/10cc syringe)

Sample gently mixed, placed on ice and timer set for **5 minutes** by: _____

500ml of 37°C RPMI with 10% AB serum added to product bag by: _____

Timer immediately set for 45 minutes by: _____ at _____

Product mixed thoroughly. 125mls from product bag dispensed into each 600ml bag containing warm media (in order of 1,2,3) and 125ml warm media dispensed into remaining product bag. The result is 4 bags containing 250mls each. Bags filled by: _____

Bags placed in incubator for remainder of the 45 minutes, mixing bags gently every 5 minutes.

Cell Processing – Post Labeling

Product removed from incubator by: _____ at _____.

250mls cold working buffer added to each bag by: _____ Each bag placed on ice for 10 minutes by: _____ (timer set)

Bags centrifuged at 300xg with no brake for 10 minutes at 4°C by: _____

Supernatant in each bag expressed by: _____ (approximately 500ml total)

200mls cold working buffer added to bag#1 by: _____. All 4 bags thoroughly mixed, rinsed and recombined into one bag by: _____

Cells centrifuged at 300xg without brake for 10 minutes at 4°C by: _____

Cell Processing – Magnetic Labeling

Supernatant expressed and cells re-suspended by: _____

Volume of cells remaining: _____ (target is smallest volume feasible)

Microbeads injected into cell mixture by: _____ (20g needle/10cc syringe)

Sample gently mixed, placed on ice and timer set for 15 minutes by: _____

After the 15 minutes incubation, bag filled with cold working buffer by: _____

Sample centrifuged at 300xg without brake for 10 minutes at 4-8°C by: _____

Supernatant removed by: _____ Volume adjusted to 100mls with cold working buffer by: _____

Sample removed for analysis by: _____

Barcode

Pre Selection Product (Sample B)

Cell count: _____ Volume: _____ Total Cells: _____ Cell Recovery: _____
 Flow sent _____ HPCA sent _____ Viability: _____

Cell Processing – Magnetic Separation

Sample collection bag prepared and weighed by: _____ Tare weight of bag: _____
 CliniMACS prepared for Enrichment 3.2 and tubing set installed by: _____

Cell Processing – Post Magnetic Separation

Weight of end product: _____ minus tare weight of cell collection bag = volume of final product: _____
 Product concentrated and reconstituted to _____ mls by _____.

Post Selection Product (Sample C – Pre-Concentration)

Cell count: _____ Volume: _____ Total Cells: _____ Flow sent on product _____,
 negative fraction _____ and buffer waste bag _____ HPCA sent _____ Product viability: _____

Post Selection Product, Post Concentration (Sample C – Post-Concentration)

Cell count: _____ Volume: _____ Total Cells: _____ Culture (21CFR.610.12) bacterial _____ fungal _____ gram stain _____ (MICRO LAB); culture _____ (STCL),
 Endotoxin _____ Flow sent on product bag _____

Supplies and Reagents

ITEM	SUPPLIER	LOT NUMBER	EXPIRATION DATE	USED
Alcohol Pads				
BacT/Alert SA				
BacT/Alert SN				
Plasma Transfer Sets (spike/spike)				
Needles (16 gauge)				
Needles (20 gauge)				
Sampling Site Coupler				
Sterile Docking Wafers				
Syringe (1ml)				
Syringe (5ml)				
Syringe (10ml)				
Syringe (20ml)				
Syringe (30ml)				
Syringe (60ml)				
Transfer Pack (150ml)				
Transfer Pack (300ml)				
Transfer Pack (600ml)				

Barcode

ITEM	SUPPLIER	LOT NUMBER	EXPIRATION DATE	USED
RPMI				
AB Serum				
25% HAS				
Gas Permeable Culture Bag				
CytoStim				
Anti-IFN γ -PE, human				
Buffer				
Tubing Set				
CMV pp65				
CCS (IFN-gamma) Microbeads				
Pall Filter				
Blood Filter				
Sterile Caps				
Plasma Transfer Set				
Luer/Spike Interconnectors				
24 Well Culture Dish				
Pre-Separation Filters				
15ml Sterile Conicals				
50ml Sterile Conicals				

Equipment

	Make	Model	Serial Number
CliniMACS			
Sterile Docker			
MiniMACS			
MACSmix			
Centrifuge			
Centrifuge			
Balance			
Tubing Sealer			
BSC			
Cell Counter			

I certify that all reagents and supplies used in processing these samples show no signs of contamination, irregularities, defects or flaws.

Date _____ Initials _____

I certify that all heat sealed tubing and all sterile dockings used in processing these samples exhibited no signs of leakage, irregularities, defects or flaws.

Date _____ Initials _____

Signature Manifest**Document Number:** STCL-PROC-037 FRM1**Revision:** 02**Title:** Miltenyi Large Scale CCB (IFN-gamma) Worksheet FRM1**STCL-PROC-037 FRM1 Miltenyi Large Scale CCB (IFN-gamma) Worksheet****Author Approval**

Name/Signature	Title	Date	Meaning/Reason
Barbara Waters-Pick (WATE02)		12 Dec 2012, 04:39:59 PM	Approved

Manager Approval

Name/Signature	Title	Date	Meaning/Reason
Barbara Waters-Pick (WATE02)		12 Dec 2012, 04:40:38 PM	Approved

Medical Director Approval

Name/Signature	Title	Date	Meaning/Reason
Joanne Kurtzberg (KURTZ001)		12 Dec 2012, 04:42:02 PM	Approved

QA Approval

Name/Signature	Title	Date	Meaning/Reason
Linda Sledge (SLEDG006)		12 Dec 2012, 05:39:27 PM	Approved

Document Release

Name/Signature	Title	Date	Meaning/Reason
Sandy Mulligan (MULLI026)		12 Dec 2012, 05:43:55 PM	Approved

Notification

Name/Signature	Title	Date	Meaning/Reason
Barbara Waters-Pick (WATE02)		12 Dec 2012, 05:43:56 PM	Email Sent
Linda Sledge (SLEDG006)		12 Dec 2012, 05:43:56 PM	Email Sent
Sharon Hartis (SH259)		12 Dec 2012, 05:43:56 PM	Email Sent
System Administrator (SYSADMIN)		12 Dec 2012, 05:43:56 PM	Email Sent