



## STEM CELL LABORATORY (STCL)



**DOCUMENT NUMBER:** STCL-PROC-042 FRM2

**DOCUMENT TITLE:**

UCB Processing Worksheet FRM2

**DOCUMENT NOTES:**

### Document Information

**Revision:** 02

**Vault:** STCL-Processing-rel

**Status:** Release

**Document Type:** SOPs

### Date Information

**Creation Date:** 05 Oct 2015

**Release Date:** 22 Oct 2015

**Effective Date:** 22 Oct 2015

**Expiration Date:**

### Control Information

**Author:** WATE02

**Owner:** WATE02

**Previous Number:** STCL-PROC-042 FRM2 Rev. 01 **Change Number:** STCL-CCR-317

## STCL-PROC-042 FRM2 UCB PROCESSING WORKSHEET

Date/Time Unit Placed on Rocker: \_\_\_\_\_

Barcode Label

Date/Time Unit Removed from Rocker: \_\_\_\_\_

**A. From Pre Processing/Collection Bag, Remove QC Sample:**

Pre Processing Sample Removed by: \_\_\_\_\_

Study ID: \_\_\_\_\_

**B. Collection Bag Calculations:**

Weight of bag = \_\_\_\_\_ mL + \_\_\_\_\_ mL CPD = \_\_\_\_\_ mL Total Volume (TV)

Study ID: \_\_\_\_\_

**C. Pre Processing Sample Analysis:**

$$\frac{\text{_____}}{\text{(WBC)}} \times 10^6 / \text{mL} \times \frac{\text{_____}}{\text{(TV)}} \text{ mL} = \frac{\text{_____}}{\text{(TNCC)}} \times 10^6 \quad \text{TNCC} = \text{_____} \times 10^8$$
**D. Processing:**OK to Process? ☐ YES ☐ NO

Study ID: \_\_\_\_\_

Total Volume / 5 = \_\_\_\_\_ mL of Hespan needed for processing

Study ID: \_\_\_\_\_

Hang Start time = \_\_\_\_\_

**Manual Processing (if applicable)** Spin @ 700 RPM for 10 min with NO brake

Let hang: \_\_\_\_\_ min Time Started: \_\_\_\_\_ Time Completed: \_\_\_\_\_ Study ID: \_\_\_\_\_

Removed: \_\_\_\_\_ mL Study ID: \_\_\_\_\_ Spin @ 1200 RPM for 20 min, then express and count

**E. Post Processing Counts/Recoveries:**

Volume Reduced/ Final Product Bag Volume = \_\_\_\_\_ mL

Study ID: \_\_\_\_\_

$$\frac{\text{_____}}{\text{(WBC)}} \times 10^6 / \text{mL} \times \frac{\text{_____}}{\text{(TV)}} \text{ mL} = \frac{\text{_____}}{\text{(TNCC)}} \times 10^6 \quad \text{TNCC} = \text{_____} \times 10^8$$

% Recovery: \_\_\_\_\_% (If recovery is below 60%, reprocess unit)

Study ID: \_\_\_\_\_

Post processing viability (%) = \_\_\_\_\_ ABO/ Rh \_\_\_\_\_

Freezer location and number of UCB plasma nuncs: \_\_\_\_\_

Freezer location and number of RBC pellet nuncs: \_\_\_\_\_

*I certify that all reagents and supplies used in processing this sample show no signs of contamination, irregularities, defects or flaws; that all heat sealed tubing and sterile welded tubing exhibit no signs of leakage, irregularities, defects, or flaws; and that all test samples related to this specific UCB unit have been identified and labeled accurately.*

\_\_\_\_\_  
Initials/Date

| Field  | Requirements  |
|--|---|
| Barcode Code Label   | Place ISBT 128 barcode for sample   |
| Date/Time Unit Placed on Rocker  | Enter date and time unit placed on rocker   |
| Date/Time Unit Removed from Rocker   | Enter date and time unit removed from rocker for processing   |
| <b>Remove QC Samples</b>   |   |
| Pre Processing Sample Removed by Study ID  | Enter study ID of person removing pre processing sample   |
| <b>Collection Bag Calculations</b>   |   |
| Weight of bag = __mL + __mL CPD = __mL<br>Total Volume (TV)  | Perform calculations for collection bag and enter your Study ID   |
| <b>Pre Processing Sample analysis</b>  |   |
| $\frac{\text{WBC}}{\text{mL}} \times \frac{\text{TV}}{\text{mL}} = \frac{\text{TNCC}}{\text{mL}} \times 10^6$ $\text{TNCC} = \text{_____} \times 10^8$ | Record WBC value from the automated count printout. Record Total Volume (TV) of collection bag and calculate Total Nucleated Cell Count (TNCC). Convert TNCC from $\times 10^6$ to $\times 10^8$ .    |
| <b>Processing OK to Process</b> <input type="checkbox"/> YES <input type="checkbox"/> NO<br>Study ID _____   | Use unit count threshold for processing per procedure STCL-PROC-042 and indicate if unit is OK to process by marking appropriate yes or no box and enter study ID                                     |
| Total Volume / 5 = _____mL of Hespan needed for processing<br>Study ID: _____  | Calculate amount of Hespan needed and enter study ID  |
| Hang time start  | Record the start time when hanging the bag.   |
| Post processing viability / ABO/Rh   | Enter viability and ABO/Rh (post processing)  |
| <b>Manual Processing (if applicable)</b>   |   |
| Let hang   | Enter the time the unit was allowed to hang   |
| Time started   | Enter time processing was started   |
| Time Completed   | Enter time that processing is completed   |
| Study ID   | Enter Study ID of person expressing cells   |
| Removed/Study ID   | Enter the amount of plasma removed from CBU and Study ID  |
| <b>Post Processing Counts/Recoveries</b>   |   |
| Volume-Reduced Final Product Bag Volume  | Enter volume from Sepax 2 S-100 printout and study ID   |
| $\frac{\text{WBC}}{\text{mL}} \times \frac{\text{TV}}{\text{mL}} = \frac{\text{TNCC}}{\text{mL}} \times 10^6$ $\text{TNCC} = \text{_____} \times 10^8$ | Record WBC value from the automated count printout. Record Total Volume (TV) of final product bag and calculate Total Nucleated Cell Count (TNCC). Convert TNCC from $\times 10^6$ to $\times 10^8$ . |
| %Recovery  | Take post TNCC and divide by pre TNCC and multiply by 100. Enter recovery % and study ID  |
| Post processing viability / ABO/Rh   | Enter viability (post processing) and ABO/Rh  |
| Freezer location and number of UCB plasma nuncs  | Enter -80 freezer storage location and number of nuncs being stored   |
| Freezer location and number of RBC pellet nuncs  | Enter -80 freezer storage location and number of nuncs being stored   |
| Certification Statement  | Initial and date to certify UCB processing  |

**Signature Manifest****Document Number:** STCL-PROC-042 FRM2**Revision:** 02**Title:** UCB Processing Worksheet FRM2

All dates and times are in Eastern Time.

**STCL-PROC-042 FRM2 UCB Processing Worksheet FRM2****Author**

| Name/Signature                  | Title | Date                     | Meaning/Reason |
|---------------------------------|-------|--------------------------|----------------|
| Barbara Waters-Pick<br>(WATE02) |       | 08 Oct 2015, 11:13:54 PM | Approved       |

**Manager**

| Name/Signature                  | Title | Date                     | Meaning/Reason |
|---------------------------------|-------|--------------------------|----------------|
| Barbara Waters-Pick<br>(WATE02) |       | 08 Oct 2015, 11:14:20 PM | Approved       |

**Medical Director**

| Name/Signature                 | Title | Date                     | Meaning/Reason |
|--------------------------------|-------|--------------------------|----------------|
| Joanne Kurtzberg<br>(KURTZ001) |       | 09 Oct 2015, 02:24:27 PM | Approved       |

**Quality**

| Name/Signature         | Title | Date                     | Meaning/Reason |
|------------------------|-------|--------------------------|----------------|
| John Carpenter (JPC27) |       | 12 Oct 2015, 09:48:56 AM | Approved       |

**Document Release**

| Name/Signature            | Title | Date                     | Meaning/Reason |
|---------------------------|-------|--------------------------|----------------|
| Sandy Mulligan (MULLI026) |       | 21 Oct 2015, 02:03:25 PM | Approved       |