



STEM CELL LABORATORY (STCL)



DOCUMENT NUMBER: STCL-SOP-058 JA3

DOCUMENT TITLE:

Entering Test Results in EPIC Beaker JA3

DOCUMENT NOTES:

Document Information

Revision: 01

Vault: STCL-Processing-rel

Status: Release

Document Type: STCL

Date Information

Creation Date: 26 Mar 2015

Release Date: 13 Apr 2015

Effective Date: 13 Apr 2015

Expiration Date:

Control Information

Author: WATE02

Owner: WATE02

Previous Number: None

Change Number: STCL-CCR-281

STCL-SOP-058 JA3

Entering Test Results in EPIC Beaker

- Log into Epic Beaker using assigned username and password.
- When the STCL dashboard is open, the tests that need to be resultated will show up the following ways:
 - On Test - Tests that have been ordered and no results have been entered.
 - Resulted - Tests that have preliminary results entered but are awaiting for results before verification
 - Pend Final - Tests that have been resultated by the STCL staff and are waiting for the interpretation to be entered by the STCL medical director.
- Double click on the test you wish to enter results.
- The test case will open, and you will need to select the “Edit” button on the left side of the screen.
- Enter the appropriate test results for that particular test. (**NOTE:** Always verify the accession number assigned for each test before entering test results, as there may be multiple tests ordered, from the same sample/product (ie. *HPC Phenotype* tests could be ordered on the same day to determine the (1) peripheral CD34 count before apheresis and to determine the (2) CD34 count from the apheresis product); these tests will have been ordered on the same day from the same patient but with different accession numbers).
- Once all of the results have been entered for that particular accession number, if an interpretation is NOT required by the medical director, that test can be verified. Select the “Verify” button located at the top center of the screen, and select “Final Verify” from the drop down window. The results entered will then appear in the window; select “Final Verify” at the bottom right of the screen.
- Make sure to double check the accuracy of all of the results BEFORE final verification.
- For any test that the medical director must verify (ie. *HPC Phenotype*, *HPCA Basic*, and *IR Panel*), enter the test results in the appropriate fields, and click “Verify”. When the drop down list appears, select “Pend Final”. When the results appear on the screen, double check for accuracy before clicking “Pend Final” at the bottom right of the screen.
- For any test that will have results entered at a later date, as is the case for bacterial cultures that take 7 days, (ie, *HPC Leukapheresis*, *HPC Reinfusion Stem Cells*, *HPC Harvest*, *Cryopreservation*), enter the test results in the appropriate fields, and select “Save” located beside the “Edit” button on the left center of screen. Once all of the final test results are available, select “Verify” then choose “Final Verify” from the drop down menu.

Signature Manifest**Document Number:** STCL-SOP-058 JA3**Revision:** 01**Title:** Entering Test Results in EPIC Beaker JA3

All dates and times are in Eastern Time.

STCL-SOP-058 JA3 Entering Test Results in EPIC Beaker**Author**

Name/Signature	Title	Date	Meaning/Reason
Barbara Waters-Pick (WATE02)		27 Mar 2015, 10:54:47 AM	Approved

Manager

Name/Signature	Title	Date	Meaning/Reason
Barbara Waters-Pick (WATE02)		27 Mar 2015, 10:54:57 AM	Approved

Medical Director

Name/Signature	Title	Date	Meaning/Reason
Joanne Kurtzberg (KURTZ001)		30 Mar 2015, 09:34:52 AM	Approved

Quality

Name/Signature	Title	Date	Meaning/Reason
John Carpenter (JPC27)		30 Mar 2015, 11:50:01 AM	Approved

Document Release

Name/Signature	Title	Date	Meaning/Reason
Sandy Mulligan (MULLI026)		02 Apr 2015, 04:32:39 PM	Approved

Quick Approval**Approve Now**

Name/Signature	Title	Date	Meaning/Reason
System Administrator (SYSADMIN)		13 Apr 2015, 09:13:08 AM	Approved